

What's Next



Centrifugation



NEW Centrifuge Bundles
Page 2

Sample Handling



Pipette / ThermoMixer Bundles
Page 3

Consumables Promotion



Buy 4, Get the 5th Box Free
Page 4

The Next Level

Innovative Bundles and Promotions for your Lab

12th August - 15th December 2019
New Zealand

Get a Bonus Spin with Eppendorf!

Now is the perfect time to start using the right rotor for your spin columns!

Centrifuge **5424 R** is the refrigerated 24-place standard microcentrifuge for temperature sensitive samples. It offers

- > Efficient cooling performance for up to 24 x 1.5/2.0 mL tubes
- > Optibowl design for quiet operation even without a lid
- > 4 rotor options, including a unique solution to prevent lids shearing from spin columns!

Bonus Kit Rotor

For a limited time, buy the 5424 R knob version with standard FA-45-24-11 rotor and get a Bonus Kit rotor free of charge!



Order No.	Bundle Description	List Price
2232001003	Bundle 5424 R, rotary knobs, FA-45-24-11 rotor, Bonus FA-45-18-Kit rotor	\$8,529.57



Eppendorf's most versatile microcentrifuge now with a Bonus plate rotor!

The 5430 R combines the best features of a refrigerated microcentrifuge (small footprint) and multipurpose centrifuge (versatility) in one instrument.

- > Max. rotor capacity: 48 x 1.5/2.0 mL, 6 x 50 mL, 2 x MTP (12 rotor options!)
- > Max. speed: 30,130 x g (17,500 rpm)
- > Soft-touch one-finger lid closure for ergonomic operation
- > FastTemp pro® function for programmable pre-cooling according to day and time

For a limited time, buy the 5430 R G, keypad version with standard FA-45-30-11 rotor & get a Bonus A-2-MTP rotor free of charge!



Order No.	Bundle Description	List Price
2230000026	Bundle 5430 R, keypad, FA-45-30-11 rotor, Bonus A-2-MTP rotor	\$11,698.32

Bonus Plate Rotor

Pump up the Volume: Bonus Offer on 5910 R Centrifuge with S-4xUniversal Rotor!

The refrigerated Centrifuge 5910 R sets the next benchmark for benchtop centrifuges.

- > Very high sample numbers per run
- > Max speed with S-4x Universal Large rotor up to 4,347 x g (S-4x400 rotor will go to 5,263 x g)
- > Excellent sample protection - Dynamic Compressor Control offers precise cooling performance at max. speed
- > FastTemp pro for automated pre-cooling programmable by time & date
- > Intuitive easy operating system and low noise level - for stress-free operation
- > Great versatility - Universal rotor, universal bucket and choice of 3 innovative universal adapters to spin a selection of tubes, bottles & plates



For a limited time, buy the 5910 R Bundle with S-4xUniversal rotor and universal adapters for 5/10/15/50 mL tubes & plates get 4 Bonus Eppendorf 1L bottles and 1L adapters free of charge!

Bonus Bottles & Adapters

Order No.	Bundle Description	List Price
2232001007	Bundle 5910 R G, S-4xUniversal rotor, universal buckets & adapters for 5/15/50 mL conical tubes & plates, Bonus Eppendorf 1L bottles (4x) and 1L adapters	\$23,616.00

Need New Pipettes?



Offer 1:		List price
Save over 30 %* when you purchase the Research plus 6-Pack Pipette and Carousel bundle		
3120000887	Research plus 6-Pack with Pipette Carousel: Includes 6 adjustable-volume pipettes (Research plus 0.1-2.5uL, 0.5-10uL, 2-20uL, 10-100uL, 20-200uL, 100-1000uL) and a Pipette Carousel 2. One box of pipette tips is included with each pipette.	\$2,150.00
Offer 2:		List price
Choose any Research plus 3-Pack value bundle & receive a free Carousel		
2232000091	Research plus 3-Pack, Option 1: Includes Research plus 0.5-10uL, Research plus 10-100uL, Research plus 100-1000uL. One box of pipette tips is included with each pipette. Includes FREE CAROUSEL	\$1,112.37
2232000092	Research plus 3-Pack, Option 2: Includes Research plus 2-20uL, Research plus 20-200uL, Research plus 100-1000uL. One box of pipette tips is included with each pipette. Includes FREE CAROUSEL	\$1,112.37
2232000093	Research plus 3-Pack, Option 3: Includes Research plus 100-1000uL, Research plus 0.5-5mL, Research plus 1-10mL. One box of pipette tips is included with each pipette. Includes FREE CAROUSEL	\$1,112.37

* Based on individual pipette RRP



Good to know – Pipetting facts and much more:
www.eppendorf.com/handling-solutions

The Value of Mixing

ThermoMixer® C and SmartBlock™

The Eppendorf ThermoMixer C offers reproducibility for temperature-controlled mixing protocols. The individually sensor-controlled Eppendorf SmartBlock provides temperature accuracy and block homogeneity for maximum reproducibility. 2D Mix-Control for controlled, circular mixing and anti-spill technology offer further assurance.



ThermoTop® and condens.protect® technology

The Eppendorf ThermoTop with special condens.protect® technology safeguards the homogeneous reaction conditions inside your tube - the system reliably prevents the formation of condensate on the tube lid and tube wall.



condens.protect®



Eppendorf ThermoMixer Bundle

Bonus

Purchase a ThermoMixer C Bundle and receive another SmartBlock or ThermoTop FREE*

2232000089	ThermoMixer Bundle with 1.5 mL Block*	\$5,483.14
2232000090	ThermoMixer Plate Bundle with ThermoTop^	\$6,696.92

* Free SmartBlock must be equal or lesser value than the SmartBlock already included in the bundle and can include any of the following: 0.5 mL / 1.5 mL / 2.0 mL / 5.0 mL / 15 mL / 50 mL / cryo blocks or ThermoTop. ^ Any free SmartBlock can be selected

Your benefits

- > Eppendorf Mix-Control ensures homogeneity throughout your sample
- > Constant reaction conditions for optimal enzymatic performance within your sample
- > Wireless connection with automatic recognition and operation for user-friendly handling
- > SmartBlock QuickRelease for easy block exchange



Special Offer

For ordering information visit:
www.eppendorf.com/consumables
 or scan the QR code



Buy 4, Get the 5th Box Free - Consumables Promotion

See it for yourself. For a limited time, take advantage of our Buy 4, Get the 5th Box free offer on a selection of consumables*.

Includes:

- > Eppendorf Tubes (includes NEW 25 mL Tube)
- > Eppendorf epT.I.P.S.[®] and Dualfilter epT.I.P.S.[®]
- > Eppendorf Combitips Advanced[®]
- > Eppendorf Plates
- > Eppendorf PCR Consumables
- > Eppendorf Serological Pipets
- > Eppendorf Cell Culture Consumables
- > Eppendorf Cell Imaging Consumables


* Free box of consumables will be supplied to you at the time of quote and time of order.

Size Matters: NEW 25 mL Tube!

Bridge the gap between 15 and 50 mL tubes with the NEW Conical Tube 25 mL

Your benefits

- > Single-handed operation with unique, reliable Snap-Tec[™] cap
- > Tight seal up to 100°C
- > Minimised contamination risk due to wider opening
- > Improved sample recovery through lower height
- > Less storage space required; less waste generated
- > A comprehensive system solution with all necessary accessories



Find our "General Quality Certificate for epT.I.P.S.[®], Tubes[®], Plates[®] and UVette[®]" here:



Not all consumables are the same

Demanding applications require careful selection of consumables which will not adversely affect the quality of analytical results. Eppendorf is committed to developing laboratory consumables to make every day routines faster, easier and more reliable. Are you aware of the many substances that leach from plastic consumables?

The most critical leaching additives are slip agents, biocides and plasticisers

False positives, wrong interpretation of biological results, waste of time and costs and unreliable measurement are some of the outcomes.

Unfortunately, leachables are easily detected in UV- absorbent measurements, interfering with spectroscopic analyses such as those employed during DNA quantifications. Eppendorf is committed to manufacturing only the highest quality of plastic consumables without the use of slip agents, plasticisers or biocides, giving you peace of mind about your results.



Authorised Distributors in New Zealand

Lab Supply Ltd.

Phone: 0800 522 787

E-mail: info@labsupply.co.nz

Website: www.labsupply.co.nz

Terms & Conditions: Offers are only valid for customers residing in New Zealand for items placed on one purchase order from August 12th - December 15th, 2019 and are available while supplies last. Promotions cannot be combined with any other offers unless otherwise stated. Free consumables must be equal or lesser value than the lowest price item purchased. All prices are in NZD and do not include packing and freight charges. For a full list of applicable products visit www.eppendorf.com/promotions. For all terms and conditions please visit www.eppendorf.com

Eppendorf South Pacific Pty. Ltd. • L1 • 97 Waterloo Road • Macquarie Park • NSW 2113 • Phone +61 2 9889 5000 • Fax No. +61 2 9889 5111 • E-mail info@eppendorf.com.au

Eppendorf[®], the Eppendorf logo, Eppendorf epT.I.P.S.[®], Eppendorf Tubes[®], Eppendorf Plates[®], Eppendorf Combitips advanced[®], Eppendorf Kit rotor[®], Eppendorf QuickLock[®], condens.protect[®], Eppendorf ThermoTop[®], Eppendorf ThermoMixer[®], and UVette[®] are registered trademarks of Eppendorf AG, Germany. Eppendorf SmartBlock[™] is a trademark of Eppendorf AG, Germany. All rights reserved, including graphics and images. Copyright © 2019 by Eppendorf AG.