


LAB
SUPPLY



IKA



**IKA – Now Exclusive to Lab Supply
in New Zealand from January 2026**

New Year Kick Off – Introductory Offer



IKA - Now Exclusive to Lab Supply in New Zealand from January 2026



IKA's product design, engineering excellence, and commitment to evolving with their customers make them a trusted partner in advancing science and industry.

Explore IKA's range of innovative laboratory solutions, designed to work perfectly and built to last.

www.ika.com

Specialising in heaters, stirrers, homogenisers, shakers, incubators, rotary evaporators, and more, IKA's portfolio is built to support a wide variety of laboratory workflows. From routine sample preparation to advanced research applications, IKA instruments are designed to deliver accuracy, repeatability, and ease of use.

What began in Cologne over 100 years ago is now an internationally recognised brand at the forefront of laboratory technology. IKA's product design, engineering excellence, and commitment to evolving with their customers make them a trusted partner in advancing science and industry.

Explore IKA's range of innovative laboratory solutions, designed to work perfectly and built to last.

labsupply.co.nz

Introductory Specials

PROMOTIONS – IN STOCK NOW

IKA



VORTEX 3
ONLY \$640



RCT BASIC
ONLY \$1,299



MATRIX Delta Plus
ONLY \$6,999



MS 3 DIGITAL
ONLY \$1,290



ROTAVISC me-vi Complete
ONLY \$6,990

TERMS & CONDITIONS

Prices for all products exclude GST. Promotional pricing applies only to products included in the IKA Introductory promotion and is not valid for all product lines. To redeem promotional pricing, customers must quote the relevant promotional code when ordering with our Customer Service team. Please note: promo codes cannot be applied online via checkout. This promotion is valid until 31 March 2026, or while stocks last. No rainchecks. The promotion cannot be used in conjunction with any other offers or discounts. Lab Supply reserves the right to amend or withdraw this promotion at any time without prior notice.

HEATER STIRRERS:

IKA

**RCT BASIC****ONLY \$1,299**

Code: IKA0025006458



RCT basic is suitable for stirring tasks up to 20 l (H₂O) and reaches a hotplate temperature of up to 310°C. With the connection option for an external temperature sensor (PT 1000.60 included in delivery), the temperature can be measured and controlled directly in the reaction medium.

HEATER STIRRERS:

IKA

C-MAG HS 7 digital**ONLY \$1,090**

Code: IKA0003581205



Magnetic stirrer with heating and ceramic heating plate which offers excellent chemical resistance.

- Powerful motor for stirring quantities of up to 10 l (H₂O)
- Simultaneous digital display of target and actual temperatures via LCD display
- Directly connection for a PT 1000 temperature sensor enable a precise temperature control (included in delivery)
- Control accuracy in medium +/- 0,5 K (in combination with PT 1000)
- Display of actual temperature value in medium with 0,1 K resolution when using PT 1000 temperature sensor

STIRRERS:

IKA

IKA MY PLATE

ONLY \$690

Code: IKA0025005156

**Personal interchangeable frame**

Design the set-up plate of the IKA myPlate by yourself. You receive the device with a decorative image template, that can be exchanged easily through the convenient interchangeable frame. Some people have family photos on their desk - you convert one of the best magnetic stirrers from your day-to-day lab work into a canvas for your favourite pictures.

- High chemical resistance
- Electronically controlled motor
- Digital speed indicator
- Speed range between 100 and 2500 rpm

SHAKER VORTEXERS

IKA

VORTEX 3

ONLY \$640

Code: IKA0003340005

**Shaking movement:**
Orbital**Orbital diameter:**
4 mm**Infinitely adjustable speed
range**
500 - 2.500 rpm

Vortex shaker suitable for short-time operation (touch function), activated by pressing shaker attachment or continuous operation.

- Different applications thanks to 3 interchangeable attachments and 7 inserts (e.g. tubes, microtiter plates, Erlenmeyer flasks 250 ml etc.), please order separately
- Attachments securely click onto appliance in any position
- Special strap (VG 3.36, page 51) ensures easy handling of round/Erlenmeyer flasks
- Sturdy metal zinc die cast casing
- Short-time operation activated by pressing attachment (touch function)
- Eccentric with ball bearings

SHAKER VORTEXERS

IKA

MS 3 BASIC

ONLY \$990

Code: IKA0003617000



Compact, universal small shaker suitable for shaking tasks with all small vessels, PCR & microtiter plates

- Wide range of attachments incl. PCR PLATE or Microplate
- Attachment detection
- Continuous or touch operation (with standard attachment)
- Two operating modes

Scope of delivery;

- MS 3 basic
- MS 3.1 Standard attachment
- MS 3.3 Universal attachment
- MS 1.21 One-hand insert

SHAKER VORTEXERS

IKA

MS 3 DIGITAL

ONLY \$1,290

Code: IKA0003319000



Compact, universal small shaker suitable for shaking tasks with all small vessels, PCR and microtiter plates

- Wide range of attachments incl. PCR PLATE or Microplate
- Attachment detection
- Timer with countdown function
- Continuous or touch operation (with standard attachment)
- Two operating modes:

Scope of delivery:

- MS 3 digital
- MS 3.1 Standard attachment
- MS 3.3 Universal attachment
- MS 3.4 Microtiter attachment
- MS 1.21 One-hand insert
- MS 1.32 Test tube insert

THERMAL SHAKERS

IKA


IKA MATRIX Orbital Delta Plus S005

ONLY \$6,999

Code: IKA0020100521



THERMAL COVER

**WITH FREE ATTACHMENT:**
you choose (more options available ie 15ml)

VISCOMETER

IKA

ROTAVISC me-vi Complete**ONLY \$6,990**Code: IKA0010013827**Scope of delivery:**

- ROTAVISC me-vi control S000
- ROTASTAND
- SP set-2
- PT 100.8
- Guard rail

The ROTAVISC series determines the viscosity of liquids in all areas of application ranging from the laboratory to quality control. The four devices measure in different viscosity ranges. Regardless of a simple or demanding viscosity measurement - the ROTAVISC delivers rapid and accurate results. The scope of delivery includes a standard spindle set (SP7- SP12), a protective bracket, temperature sensor and ROTASTAND stand.

Easy to operate

The generous 4.3" TFT display allows for an intuitive menu guidance. A digital spirit level also supports the correct set-up of the start-up condition. The simple ramp function simplifies repetitive tasks.

Highest measurement accuracy

The level of measurement accuracy of the ROTAVISC rotational viscometer for both Newtonian and non-Newtonian fluids is +/- 1 % of the measuring range. The reproducibility is +/- 0.2 %.

Stepless speed

ROTAVISC is a viscometer, which offers stepless speed adjustment.

IKA ROTAVISC me-vi Complete

Viscosity measuring range: 100 - 40.000.000 mPas

SHOP IKA's full range at LAB SUPPLY



Rotary Evaporators

Sets new standards for safety, efficiency and ease of use



Overhead Stirrers

Effortlessly stir up to 200 litres at any viscosity



Bomb Calorimeter

Sets new standards for safety, efficiency and ease of use



Homogenisers

Know no limits when it comes to professional particle reduction



Mill Range