



Flex-ibility and affordability in the perfect bundle.

See your lab's budget in a whole new light with our Reflectoquant® system bundle!

Maintain your high standards even when you are on-the-go with our Reflectoquant® system. The quantitative test strip reader can monitor your chemical parameters on the spot.

Our RQflex® 20 reflectometer enables you to take faster and traceable decisions in your production and disinfection control process for rapid quantitative results—directly on-site.

Take advantage of our amazing bundle offer to get the RQflex® 20 reflectometer **PLUS** a selection of test kits at an unbeatable price!


So don't reflect on this for too long - get moving before this offer expires! Find out more on the next page, or contact your account manager today.

To learn more about the Reflectoquant® system visit:

SigmaAldrich.com/Reflectoquant

Available from:

LABSUPPLY 



Reflectoquant® Bundle Information*

Get **1** each of:

Product	Content	Cat. No.
Reflectometer RQflex® 20 Reflectoquant®	Includes test strip adapter and recalibration set, double optical system (option for evaluation of two reaction zones), memory for five methods, memory slots for 50 results (with date, time, parameter, and result), batch-specific calibration function (barcode technology), battery operation with four 1.5 V batteries, quick guide for reflectometer and tests	1.17246.0001
RQcheck set for RQflex 20 Reflectoquant®		1.17247.0001

AND Any **3** of:

Reflectoquant® tests	Measuring range	Number of tests	Cat. No.
Recalibration set for RQflex® 20			1.16954.0001
Nitrate Test	3 – 90 mg/L NO ₃	50	1.16995.0001
Nitrate Test	5 – 225 mg/L NO ₃	50	1.16971.0001
Nitrite Test	0.5 – 25 mg/L NO ₂	50	1.16973.0001
Ascorbic Acid Test	25 – 450 mg/L ascorbic acid	50	1.16981.0001
Calcium Test	5 – 125 mg/L Ca	50	1.16125.0001

OR Any **2** of:

Reflectoquant® tests	Measuring range	Number of tests	Cat. No.
Total Sugar Test (glucose and fructose)	65 – 650 mg/L total sugar	50	1.16136.0001
Nitrite Test	0.03 – 1.00 g/L NO ₂	50	1.16732.0001
Ammonium Test	20 – 180 mg/L NH ₄	50	1.16977.0001
Hydroxymethylfurfural Test	1.0 – 60.0 mg/L HMF	50	1.17952.0001

OR **1** of:

Reflectoquant® tests	Measuring range	Number of tests	Cat. No.
Sucrose Test	0.25 – 2.50 g/L	50	1.16141.0001

For a full list of included Reflectoquant® tests, get your copy of the WFA Catalogue today:

SigmaAldrich.com/WFA-catalog

Contact your local distributor for sales support:

Lab Supply Ltd
Ph: 0800 522 787
Email: info@labsupply.co.nz
Web: www.labsupply.co.nz



**Bundle
Price:**
\$1190_{ex-GST}
T's & C's Apply*

*Terms & Conditions:

Offer Valid until 31st December 2023 for Australia and New Zealand customers only.

Bundle includes one 1172460001 RQflex® 20 unit plus one 1172470001 RQ Check Set and a combination of Reflectoquant® consumables including but not limited to those shown in the tables above. Acceptable kits and quantities will vary based on individual list price.

Bundle price total \$1190 NZD ex-GST.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, Supelco, Spectroquant and the vibrant M are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

The life science business
of Merck operates as
MilliporeSigma in the
U.S. and Canada.

Supelco®
Analytical Products